



Patient information: pharmacy

Compliance aids step by step guides-Alcon Eyot®

The Alcon Eyot® is for eye drop bottles made by the manufacturer Alcon, fitting their specific shaped bottles- Travatan® Duotrav® and Azarga®.

It is designed to help with the stable positioning of the bottle over your eye, so you can accurately squeeze a drop into your eye.

Instructions for use:

- After washing your hands, insert your eye drop bottle into the grooves by the mirror (found inside the aid), with the nozzle pointing up towards the mirror.
- 2. Lie down or sit down comfortably, then twist off the bottle cap.
- 3. Tilt your head back and place the aid over your eye so that it pulls your bottom lid down to form a pocket.
- 4. Look in the mirror. When you can see your eye, the bottle should be in the right position.
- 5. Look through the pin hole at top, and when ready, squeeze the sides to release a drop into your eye.
- 6. Unless your clinician has advised otherwise, close your eye and gently press the inside corner for 30 seconds to help prevent the drop entering your mouth.
- 7. When finished, replace the cap of your eye drop bottle.

Author: Sarah Thomas, lead satellite pharmacist and Fiona Chiu, associate chief pharmacist

Revision number: 1 Approval date: April 2018

This information can be made available in alternative formats, such as easy read or large print on request. Please call PALS: 020 7566 2324 / 020 7566 2325

